



Thumb and Wrist/Thumb Spica

HCPCS Code L3807



- Lightweight splint constructed of perforated foam
- Thumb spica supports thumb while allowing IP joint to remain free
- Wrist and thumb spica supports the wrist and IP joint of the thumb
- Soft inner lining wicks away perspiration for comfort and hygiene
- Stretch loop straps designed for a precise anatomical fit allowing hand and finger function
- May be trimmed for a custom fit

INDICATIONS: De Quervain's Syndrome; CMC arthritis; Post-thumb fractures; Tendon injury/repair

Thumb

Size	Measurement	Right	Left
Small/Medium	2½" - 3½"	340S/MR	340S/ML
Medium/Large	3½" - 4¼"	340M/LR	340M/LL

Thumb/Wrist

Size	Measurement	Right	Left
Small/Medium	2½" - 3½"	345S/MR	345S/ML
Medium/Large	3½" - 4¼"	345M/LR	345M/LL