

Wrist-Hand PlastiCast

EMAIL A FRIEND 



PRODUCT HIGHLIGHTS

- Immobilizes wrist while allowing full finger/thumb ROM
- Semirigid polyethylene shell with a 1/8" thermo-bonded closed-cell polyethylene foam interior for comfort and easy application

MORE ABOUT THIS PRODUCT

Wrist-Hand PlastiCast immobilizes wrist while allowing full finger/thumb ROM. Semi-rigid polyethylene shell with a 1/8" thermo-bonded closed-cell polyethylene foam interior for comfort and easy application. Hook and loop closures. Ventilated. Hypoallergenic. Radiolucent. Water resistant. Modify with scissors and heat gun. Comes with one pre-sewn, cotton stockinette liner.

Indications: Post-cast fracture immobilization, navicular fractures, nondisplaced distal forearm fractures.

Suggested code: L3908.

Sizing: MP diameter, Wrist Circumference, Forearm Circumference

Sizing	MP Dia.	Wrist Circ.	Forearm Circ.
Small	2-3/4" to 3-1/4"	7" to 8"	8" to 10"
Medium	3-1/4" to 3-3/4"	8" to 9"	10" to 12"
Large	3-3/4" to 4"	9" to 10"	12" to 14"